

## DP-6900

### Digital Ultrasonic Diagnostic Imaging System



#### Features:

- ▶ Tissue Harmonic Imaging (THI), Tissue Specific Imaging (TSI)
- ▶ Multipurpose----9 transducers available with up to 6 frequencies
- ▶ Comprehensive exam modes and measurement packages (OB/GYN, abdomen, small parts, cardiology, urology, orthopedics, vascular, breast, prostate, MSK, nerve, intraoperative), supporting user-defined functions
- ▶ iStation™-----Patient information management system
- ▶ User-defined soft keys, one key quick image storage to local disk or USB drive
- ▶ Multi-language display (English, French, German, Italian, Russian, Portuguese, Spanish and Chinese)

#### Standard Configuration:

- ▶ DP-6900 main unit
- ▶ 10" non-interlaced monitor
- ▶ THI (Tissue Harmonic Imaging)
- ▶ iStation™ Patient Information Management System
- ▶ Two USB ports
- ▶ 400M standard image-storage capacity
- ▶ Measurement & calculation software packages
- ▶ One transducer connector

**Optional:**

- ▶ Electronic convex array transducer: 35C50EA (2.0/3.5/5.0/6.0/H4.6/H6.0MHz)
- ▶ Electronic linear array transducer: 75L38EA (5.0/7.5/8.5/10MHz)
- ▶ Electronic linear array transducer: 75L60EA (5.0/7.5/8.5/10MHz)
- ▶ Electronic linear array transducer: 75L53EA (5.0/7.5/8.5/10MHz)
- ▶ Electronic micro-convex array transducer: 35C20EA (2.0/2.5/3.5/5.0MHz)
- ▶ Electronic micro-convex array transducer: 65C15EA (5.0/6.5/7.5/8.5MHz)
- ▶ Electronic endocavity transducer: 65EC10EA (5.0/6.5/7.5/8.5MHz)
- ▶ Electronic intrarectal linear array transducer: 65EL60EA (5.0/6.5/7.5/8.5MHz)
- ▶ Electronic intraoperative T-shape transducer: 75LT38EA (5.0/7.5/8.5/10MHz)
- ▶ Second transducer connector
- ▶ 80G massive image-storage module
- ▶ Needle-guided brackets
- ▶ DICOM 3.0
- ▶ Mobile trolley
- ▶ Footswitch
- ▶ Hand carried bag
- ▶ Transportation & operation bag
- ▶ Rechargeable battery



**Multi-frequency Transducers:**

---



Applications: Abdomen, GYN, OB, Urology, Orthopedics, Breast, Musculoskeletal, Endovaginal, Endorectal, Small Parts, Pediatric, Cardiac



**Gallery:**

carotid

cervical plexus

Head of Fetus

IUD

Kidney Stone

Pancreas

Liver

Thyroid Tumor

Tendon

Uterus

Thyroid

Please notice: The quality of images presented on the World Wide Web varies greatly depending on the computer system, monitor, and web browser used to visualize the images. For that reason, Mindray does not recommend that you use the images as shown here to judge the system image quality. Write to Mindray to receive hard copy images for this purpose.