



CAPNOCHECK® II

The BCI, Inc. Capnocheck®II Hand-Held Capnograph/Oximeter is a full-featured monitor providing ETCO₂, Inspired CO₂, Respiration Rate, SpO₂, and Heart Rate measurements.

This highly portable unit is ideal for use in and out of the hospital environment. A back-light feature on the display enables easy viewing in difficult lighting. Battery operated with an optional battery bypass offers flexibility for stationary and transport applications.

- ▲ Waveforms, numeric values and on-screen trending provided on an easy-to-read EL back-lit LCD display
- ▲ Adjustable alarm limits as well as visual and auditory alarm signals
- ▲ Sidestream technology for applications on intubated and non-intubated patients
- ▲ Lightweight for easy transport
- ▲ Infrared link to an external printer for trouble-free downloading

